
## ENGLISH

The Novalis conical collimators are available in various different sizes.

**Warning: Do not mount a Novalis conical collimator if the patient is directly beneath the linac gantry.**

### Exchanging the Novalis Conical Collimator

To exchange a Novalis conical collimator e.g. for another Novalis conical collimator with a different aperture:

1. Remove the fixation ring from the m3 Novalis conical collimator mount as shown in figure 1 above.

**Warning: Make sure that the security pin is correctly snapped in before removing the fixation ring. Only in this position can it be guaranteed that the Novalis conical collimator remains securely attached, so that the fixation ring can be removed without the Novalis conical collimator falling off.**

2. Firmly hold the Novalis conical collimator and press the security pin of the Novalis conical collimator towards the cone as shown in figure 2 above. This opens the locking mechanism.
3. Rotate the Novalis conical collimator counterclockwise to remove it.
4. To insert a new Novalis conical collimator, align the collimator to be mounted so that its three protrusions match the three spacing inserts within the m3 collimator mount, as shown in figure 3 above.
5. Insert the Novalis conical collimator and turn it clockwise until the security pin snaps in.
6. Replace the fixation ring and fasten it securely, as shown in figure 4 above.

**Warning: If the screws of the m3 collimator mount are opened or the m3 collimator mount is removed, the m3 collimator mount must be realigned using the Winston-Lutz test. For further information, contact BrainLAB support.**

**Warning: Make sure that the Novalis conical collimator is securely fastened before beginning patient treatment.**

**Warning: Operate the LINAC with photons only when using Novalis conical collimators.**

### Jaw Positions

It is important to verify by inspection that the area outside the shielded central part of the Novalis conical collimator is completely covered by the jaws. The models and range of Novalis conical collimator apertures differ significantly at each customer site.

For Novalis conical collimators, to avoid radiation outside the Novalis conical collimator or leakage at the edge of the field resulting from potential set-up inaccuracies (mechanical Novalis conical collimator alignment, jaw positioning, and other tolerances), BrainLAB recommends using the jaw settings listed in the table below. If properly implemented, these settings will aid in keeping the jaw shaped field size tight around the largest Novalis conical collimator.

Regardless of the jaw-shaped field size in use, BrainLAB recommends using a suitable film test to verify and document the complete shielding outside the intended center beam through the Novalis conical collimator.


Please keep in mind that the jaw-shaped field must be symmetric to the central beam axis for cone treatments, and the identical jaw position must be used for all Novalis conical collimators. If you intend to modify the jaw positions after the initial beam data acquisition you should determine to which degree the adjustment of jaw positions requires re-measurement of certain beam profiles.

Novalis (non-Novalis Tx)	
Diameter of <u>largest</u> conical collimator in use at site [mm]	Recommended size of jaw-shaped field for all conical collimators [mm <sup>2</sup> ]
Up to 7,5	14x14
8	14x14
10	16x16
12,5	18x18
15	20x20
17,5	24x24
20	26x26
22,5	28x28
25	30x30
27,5	34x34
30	36x36
32,5	38x38
35	40x40
37,5	40x40
<b>Warning: The use of a 50x50 jaw-shaped field for Novalis (non-Novalis Tx) may result in minor dose leakage at the corners of the jaw area.</b>	

### Circular Field Size

The engraved cipher on the Novalis conical collimator specifies the circular field size  $\phi$  (mm) in the isocenter.

### Cross-section of the Novalis Conical Collimator

**Warning: For Novalis conical collimators the LINAC jaws have to shield the radiation outside the cone material, but also may not be too narrow to shield the intended inner radiation field defined by the Novalis conical collimator. The LINAC jaw positions used during treatment must be identical to the LINAC jaw positions used during beam measurements.**

**Warning: To avoid unintended leakage radiation using Novalis conical collimators the LINAC jaw size must be smaller or equal 40x40 mm<sup>2</sup> for Novalis (non-Novalis Tx).**


**Warning: Ensure that the LINAC jaws do not overlap the circular Novalis conical collimator opening.**

**Warning: Be aware that the positioning accuracy of the LINAC jaws can vary depending on the type of LINAC used. The LINAC jaw positions used during treatment must be identical to the LINAC jaw positions used during beam measurements.**

**Warning: Verify for every treatment field that the correct, intended Novalis conical collimator is actually mounted to the LINAC, especially if a treatment plan contains treatment fields with different Novalis conical collimator sizes.**

**Warning: Ensure that any treatment accessories, such as m3 collimator mount and Novalis conical collimator of the size which is required by the treatment plan, are installed properly in the beam path before delivery. Also refer to the documentation provided for the BrainLAB Treatment Planning System.**

For further information please see BrainLAB's user manual or consult BrainLAB customer support.

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