

ADVANCING  
CANCER  
TREATMENT

# MLC MODELING IN RAYPHYSICS

---

RayStation 11B



---


# PROJECTIONS

- All distances for the different parameters are the distances projected on the isocenter plane.

EDUCATIONAL MATERIAL

# AUTO-MODELING

- The auto-modeling adjusts parameters to make computations match data better.
- Relevant generalizations outside the data in beam commissioning must be verified thoroughly before the model can be used clinically.



---


## TWO OPTIONS IN MLC MODELING

- Two possible workflows for MLC modeling, either importing arbitrary fields or work in the Beam 3D modeling module and external.
- From Raystation 10A, it is possible to import arbitrary fields to the Beam commissioning module to model the MLC parameters.
- MLC modeling should always be performed with any of the suggested workflow before commissioning for clinical use



# ARBITRARY FIELDS IN BEAM COMMISSIONING

- Import of arbitrary fields allows for modeling of the MLC parameters directly in the beam commissioning module.
- Possible to add any photon fields with couch and gantry angle 0°
  - No wedge, blocks or cones supported
  - No depth dose curves needed
  - Supports 3DCRT, SMLC, DMLC fields
- The list of fields for arbitrary fields is stored in the patient database and can only be accessed with scripting.
- Not possible to import dose curves for arbitrary fields in \*.xmcdat format.



# ARBITRARY FIELDS IN BEAM COMMISSIONING

- All dose curves belonging to the same arbitrary field must have the same SSD, but can have a different SSD than the rectangular fields.
- Several depths and curve types can be imported for one field.
- Each dose curve needs a scale factor that converts the measured dose to cGy.
- The arbitrary fields are not included in machine backups.
- Important that fields contain all the relevant information on the MLC property that is modeled and do not depend on several parameters for proper modeling, such as picket fence fields.

## Scale factors

*Scale factors for arbitrary fields. Dose values in the imported dose curve files will be scaled with this factor. Enter a factor such that dose values in the curve file multiplied with this factor correspond to the delivered dose in cGy. Output factor depth is not used for arbitrary fields.*

# ADMINISTRATION OF THE ARBITRARY FIELD LIST

- Can be done either through Raystation or Beam 3D modeling.

- In the command console:

- **Add a beam from a plan:**

patient\_db.CopyBeamToBeamCommissioningField(BeamToCopy= beam, FieldName='name', Description='description of field')

- **Remove a beam:**

patient\_db.RemoveBeamCommissioningField(FieldName = 'name')


- **List all beams in field list:**

patient\_db.ListAllBeamCommissioningFieldNames

- **More Information in**

**RayStation 11B RayPhysics Manual → AppendixC**

- After being added into the list, measurements can be imported to the beam model in RayPhysics



```
RayStation Command Console - "C:\Program Files (x86)\IronPython 2.7.1\ipy64.exe" -X:Tab...
File Edit View Help


IronPython 2.7.1 (2.7.0.40) on .NET 4.0.30319.42000
Type "help", "copyright", "credits" or "license" for more information.
Connecting to RayStation. (Session id = 6272_12 )
Variable 'patient' is set to patient with name: "DKFZ_LargePhantom_Pyramids^RayC
ompare - Do not edit!"
Variable 'case' is set to case with name: "CASE 1"
Variable 'plan' is set to plan with name: "InvPyramid SSD90 6MV"
Variable 'beam_set' is set to beam set with name: "test plan1"
Variable 'examination' is set to examination with name: "CT 1"
Variable 'clinic_db' is set to the current clinic database
Variable 'machine_db' is set to the current machine database
Variable 'patient_db' is set to the current patient database
Variable 'ui' is set to the current ui
Imported statetree
>>> b= beam_set.Beams[0]
>>> patient_db.CopyBeamToBeamCommissioningField(BeamToCopy=b, FieldName= 'Invers
e Pyramid', Description= 'Test')
<ScriptObject id=b0c2f675-0942-4312-a8c2-572be69f91e9, 'Inverse Pyramid'>
>>>
```

# TRANSMISSION PARAMETER

---



# TRANSMISSION





Open: transmission 1

Under Y-jaw: transmission 0


Under MLC: transmission  $t_{MLC}$

Under MLC + X-jaw:  $t_{MLC} * t_{X-jaw}$

Leaf tip or T&G: transmission  $\sqrt{t_{MLC}}$



# MLC TRANSMISSION PARAMETER




- Transmission measured on both leaf bank at 1.5 cm depth


# EXAMPLES ON TEST FIELDS FOR MLC MODELING

---

# OFFSET, GAIN AND CURVATURE



Model without correction      Reality (example)




$x_{effective} = x - offset + gain \cdot x - curvature \cdot x^2$       left MLC bank

$x_{effective} = x + offset + gain \cdot x + curvature \cdot x^2$       right MLC bank

	Positive	Negative	Effect	Range
Offset	Shifts both sides outwards	Shifts both sides inwards	All field sizes are shifted equally much	-0.0X to 0.0X
Gain			The shifts increases linearly with field size	-0.00X to 0.00X
Curvature			The shifts increases quadratic with field size	-0.000X to 0.000X


# OFFSET, GAIN AND CURVATURE



- Profiles of square MLC-collimated fields taken at 5, 10 and 20 cm depth
- Field dimensions 2x2, 3x3, 4x4, 5x5, 10x10, 15x15, 20x20 and 30x30 cm<sup>2</sup>
- Sources modeled with jaws only collimated fields



# EXAMPLES RESULTS – OFFSET, GAIN AND CURVATURE



- Optimal values: Offset: 0.04 cm, Gain: 0.0 cm<sup>-1</sup>, Curvature: 0.0cm<sup>-2</sup>

# LEAF TIP WIDTH

- Rounded lead tip model

- Region in cm (at isocenter plane), cut from the leaf tip for which the transmission equals the square root of the full leaf transmission






Fig. 11. Radiograph showing the high dose line that appears when a film is double exposed with the leaves abutted at the field midline. Notice also the diamond shaped pattern that results from the beveling at the front face of the leaves.



## EXAMPLE: LEAF TIP WIDTH



- Create an abutment using the leaf tip
- Profiles taken at 10 cm depth




# EXAMPLE RESULT – LEAF TIP WIDTH - BEAM 3D MODELING

Optimal value: leaf tip = 0.3 cm




# TOUNGE AND GROOVE

- Tounge and Groove model
  - Region in cm (at isocenter plane), extended on the leaf side for which the transmission equals the square root of the full leaf transmission









# EXAMPLE: TOUNGE AND GROOVE REGION



- Create an abutment using the leaf side
- Profiles taken at 10 cm depth

# EXAMPLE RESULT – TOUNGE AND GROOVE BEAM 3D MODELING

Optimal value: TG = 0.05 cm



# RAYSTATION DATA FROM GLENN ET AL

		Tongue and Groove		Leaf tip width		MLC transmission		MLC Offset		MLC gain		MLC Curvature		Sample size
		Mean	Range	Mean	Range	Mean	Range	Mean	Range	Mean	Range	Mean	Range	
Energy	Type													
6	Elekta Agilty	0.065	0-0.11	0.293	0.1-1	0.0043	0-0.016	0.0321	-0.02-0.12	0.0187	0-0.075	0.0175	-0.0006-0.08	40
10	Elekta Agilty	0.0792	0-0.5	0.285	0.01-0.7	0.0038	0-0.01	0	-0.02-0.02	0.0017	0-0.0047	0.0001	0-0.0004	23
15	Elekta Agilty	0.0569	0-0.1	0.239	0.1-0.65	0.0046	0-0.007	0.0565	-0.015-0.1	0.0429	0.003-0.075	0.0539	0-0.1	17
18	Elekta Agilty	0.06	0.05-0.1	0.086	0.01-0.2	0.0145	0.001-0.0365	0.0212	0.003-0.03	0.00322	0.0012-0.0085	-0.000066	-0.0006-0.00009	5
6 FFF	Elekta Agilty	0.053	0.05-0.08	0.1408	0.008-0.4	0.00276	0.001-0.005	0.017	0.01-0.025	-0.00166	-0.0046-0.0009	0	0	9
10 FFF	Elekta Agilty	0.0375	0-0.05	0.2025	0.01-0.5	0.005	0.001-0.0095	0	0	0	0	0.000175	0-0.00035	8
6	Elekta MLCi	0.064	0-0.15	0.22	0.2-0.3	0.0094	0.005-0.012	0.0233	0.0-0.05	0.0033	0-0.008	0.0004667	0.0002-0.0007	5
10	Elekta MLCi	0.025	0-0.05	0.2	0.2-0.2	0.0075	0.005-0.01	0.02	0.02-0.02	0	0-0	0.00035	0.00035-0.00035	2
15	Elekta MLCi	0.025	0-0.05	0.2	0.2-0.2	0.0075	0.005-0.01	0.02	0.02-0.02	0	0-0	0	0-0	5
6	Millenium 120	0.0396	0.01-0.05	0.2744	0.165-0.5	0.01723	0.007-0.025	0.03463	0-0.12	0.0028	-0.0004-0.015	0.000448	0-0.001	27
10	Millenium 120	0.0353	0.01-0.05	0.259	0.01-0.5	0.0181	0.01-0.0294	0.028	0-0.1	0.008	0-0.02	0.00016	0-0.0008	7
15	Millenium 120	0.0425	0.02-0.05	0.24125	0.165-0.3	0	0-0	0	0-0	0	0-0	0.0004	0-0.0008	4
18	Millenium 120	0.045	0.02-0.05	0.2275	0.165-0.3	0.0192	0.015-0.025	0.0456	0-0.12	0.00124	0-0.0031	0.00061	0-0.0008	6
6 FFF	Millenium 120	0.03889	0.01-0.05	0.6111	0.2-1	0.008878	0.0005-0.02	0.0115	0-0.042	-0.000375	-0.0015-0	0.00083	0-0.00166	9
6	Millenium HD-120	0.05	0.05-0.05	0.3	0.2-0.4	0.01627	0.01254-0.02	no data	no data	no data	no data	no data	no data	2
10	Millenium HD-120	0.05	0.05-0.05	0.25	0.2-0.3	0.02175	0.0141-0.0294	no data	no data	no data	no data	no data	no data	2

Glenn et al.: Reference dataset of users' photon beam modeling parameters for the Eclipse, Pinnacle and RayStation treatment planning systems. Medical Physics 47 (1), January 2020

EDUCATIONAL MATERIAL

# RAYSTATION DATA FOR VARIAN MACHINES FROM JORDI SAEZ ET AL 2020

- Proposed workflow to determine MLC values from: Jordi Saez *et al* 2020 *Phys. Med. Biol.* **65** 155006

**Table 1.** Summary of the optimal MLC configuration parameters obtained with the proposed procedure.

	MLCHD		Millennium		
	6 WFF	6 FFF	6 WFF	10 WFF	15 WFFF
$T$ (%)	1.19	1.03	1.46	1.69	1.64
$w_{TG}$ (mm)	0.44	0.43	0.44	0.44	0.44
$l_{tip}$ (mm)	1.19	1.05	1.64	1.81	1.85
$x_{off}$ (mm)	0.12	0.12	0.46	0.50	0.49
DLG (mm)	0.47	0.41	1.28	1.41	1.40
$k$ (Gy/(MU mm <sup>2</sup> ))	0.262 9	0.248 2	0.132 5	0.147 3	0.154 4

Image: Jordi Saez *et al* 2020 *Phys. Med. Biol.* **65** 155006

# ASYNCHRONOUS SWEEPING GAPS

- Users meeting 2020.
    - Presentation in community.
  - Publication can be found by link from RaySearch webpage (or Google).
  - Plans, scripts and analysis excel sheet provided upon request from RaySearch or the authors.
- 
- Workflow to find MLC parameters
  - Instead of high spatial accuracy work with average doses.
  - In x by sweeping different gaps.
  - In y by having a large detector, Farmer.

## Assessment and improvement of MLC models

Victor Hernandez  
vhernandezmasgrau@gmail.com



Jordi Saez  
jordi.saez@gmail.com



Saez J, Hernandez V, Goossens J, De Kerf G, Verellen D, A novel procedure for determining the optimal MLC configuration parameters in treatment planning systems based on measurements with a Farmer chamber, Phys. Med. Biol. 2020;(65):155006.

EDUCATIONAL MATERIAL




# ASYNCHRONOUS SWEEPING GAPS

- Sweeping gaps
  - Dose versus gap size gives dose to be added by leaf tip and offset
- Asynchronous sweeping gaps, aSG.
  - Fixed gaps, different extension of every second leaf, s.




- Farmer chamber measurement in solid water. Per beam.
- Extract computed dose using script. Per beam.




# ASYNCHRONOUS SWEEPING GAPS

- Increasing  $s$  -> increases shadowing -> decreases dose.
- RayStation
  - While  $s < w_{LT}$ , we are only exposing the corner giving no dose decrease.
  - When  $s > w_{LT}$  open and leaf regions are replaced with TG-regions. Linear dose decrease.



- Reshuffle measured dose to transmission reduction at one leafside .  $\Delta\Phi(s) = (D(s=0)-D(s))/2k$ .

(a)



(b)







FIG. 5: (a) Experimental fluence reduction for the aSG tests obtained for different sweeping gap sizes for the 6FFF beam and the HDMLC. (b) Experimental fluence reduction for the aSG test with a 20 mm gap and the corresponding fit for the MLC model in RayStation.

# FINDINGS

- Same shape for same treatment device
  - MLC, energy, linac type, FF/FFF.
- Varians – constant + linear fit is good
- Agility
  - Tilted TG + rounded LT gives dose reduction curve where it is not obvious which is the best constant + linear fit.
  - Comb/TG fields fits the region far from tip.
- Offset must be tuned to the LINAC.
  - First fix T, wLT and wTG.
  - Fine tune offset using patient plans but with ionization chambers rather than only PSQA-device.
    - Independency
    - Higher dose accuracy


Saini et al. Unlocking a closed system: dosimetric commissioning of a ring gantry linear accelerator in a multivendor environment, Rad. Onc. Phys. 22(2021) 21-34.

Koger et al., Impact of the MLC leaf-tip model in a commercial TPS: Dose calculation limitations and IROC-H phantom failures, Rad. Onc. Phys. 21 (2029) 82-88.




# MODELING ELEKTA CORNER BEHAVIOR

- Standard Elekta



- Elekta Agility



---

## MODELING ELEKTA CORNER BEHAVIOR

- Three options...

EDUCATIONAL MATERIAL



## OPTION 1 – SET MAXIMUM LEAF POSITION IN MLC TAB (RECOMMENDED)


From v10A

- Possible to set minimum and maximum leaf position for individual leaves in the MLC tab.
- Possible to properly model the Elekta corners behavior by accurately constrain the leaf position.
- This is the recommended way to model the Elekta corners

	Leaf center position [cm]	Width [cm]	Minimum tip position [cm]	Maximum tip position [cm]	
1	-19.75000	0.5000	-15.00	16.10	M
2	-19.25000	0.5000	-15.00	16.70	M
3	-18.75000	0.5000	-15.00	17.30	Sc
4	-18.25000	0.5000	-15.00	17.80	Th
5	-17.75000	0.5000	-15.00	18.30	M
6	-17.25000	0.5000	-15.00	18.80	M
7	-16.75000	0.5000	-15.00	19.20	
8	-16.25000	0.5000	-15.00	19.70	

# OPTION 2 – CONSTRAIN THE BEAM PROFILE CORRECTIONS

From v9A




- Model the BPC using diagonal measurements.
- Add a point in the Beam Profile at the relevant distance and set it to zero.
- Might be destroyed by auto modeling.
- Diagonal profile should always be included even though the leaf constraint is possible from v10A

EDUCATIONAL MATERIAL

# OPTION 3 - SCRIPTING


## Previous versions



- The script will examine all fields and indent the MLC in relevant fields.
- Dose will be recomputed and the plan quality can be evaluated.
- Drawbacks: Not possible to optimize with the MLC indented.

# COMPARISON FOR AN ELEKTA AGILITY HEAD



90cm SSD 5cm deep diagonal profile for 40x40cm  
6MV field



- Beam Profile Corrections for radius > 20 cm must be modelled in the "Beam 3D Modelling" tab comparing calculated and measured diagonal profiles
- Tuning the beam model fits better with measurements compared to the script


# MODELING VARIAN MACHINE CORNERS

- Different Varian machine models can have a different behavior at the corner
- The cutoff for Varian machines is due to the primary collimator




## ELEKTA AGILITY

- SSD = 100 cm Cut Off = 24.8 cm



## STANDARD ELEKTA

- SSD = 100 cm Cut Off = 24.7 cm







# BEAM PROFILE CORRECTIONS AND OFF-AXIS SOFTENING FOR RADIUS > 20 CM

- Beam Profile Corrections and Off-axis softening for radius > 20 cm – import diagonal profiles and model in “Beam Commissioning” tab
- **Note:** coordinates should be changing in both X and Y

Tank 0°  
Collimators 0°  
Diagonal measurements



*As of V9A  
and over!*



THANK YOU!